MODEL 1170-A

TR

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
11. Elektrik Devre Şemaları

11.1. Model 1170-A Elektrik Devre Şeması

1. Introduction	10
1.1. Intended Use	10
2. Description and Technical Specifications	11
2.1. General Specifications	11
3. Packaging, Storage and Transportation	12
4. Assembly and Intallation	12-13
4.1. Installation of Body	12
4.2. Assembly of Upper and Lower Pipes	12
4.3. Assembly of Spring Arm	12
4.4. Spring Arm Adjustment	12
4.4.1. Balance Adjustment	12
4.4.2. Height Adjustment	12
4.5. Assembly of the Light Head	13
4.6. Removing and Axing the Sterilisable Handle	13
5. Use	13-15
5.1. Preparation for Use	13
5.2. Description of the Control Panel	14
5.3. Adjusting the direction of light	14
5.4. Control keys	14
5.5. Opening / Closing the Lamp	16
5.6. Adjustment	14
5.7. Large / Small Lighting Focus	15
5.8. Endo Light	15
5.9. Activating / deactivating the sensor	15
5.10. Activating / deactivating the sensor	15
6. Maintenance / Cleaning And Disinfecting	15-16
6.1. Cleaning/Disinfecting	16
6.2. Cleaning the Transmissive Surface	16
6.3. Handle Cleaning and Sterilization	
Basic knowledgeg	16
7. Labeling Information	17
8. Disposal	17
9. Technical Specifications	17
10. Trouble-shooting & Service	18
11. Electrical Circuit Diagrams	18
11.1. Electrical Circuit Diagram for Model 1170-A	18

1. Introduction

Thank you for preferring BIÇAKCILAR brand.

Each device that bears BIÇAKCILAR brand reflects a production experience of more than 65 years. The reason behind the leading position of BIÇAKCILAR both in domestic and international markets is the company's uncompromising firm stand concerning the basic principles and consistency of the quality policy. "SAFETY" is the most emphasized principle at BIÇAKCILAR. It is of extreme importance that the contents of this user's manual be strictly followed to ensure the safe operation of the device. Otherwise, hazardous situations effecting the user, patient, or the environment may arise.

Important Warning!

1. Always read the user's manual prior to installation, assembly and operation of the instrument.
2. The customer is responsible for informing the relevant people on the content of this manual.
3. Be sure that the content of the user's manual is well understood prior to any assembly and operation.
4. The user's manual should be kept ready for the users.
5. In order to prevent electric shock, surgical lights must be connected to the main power source only by protective grounding.
6. The electricity system at the site where the instrument is to be installed should comply with TS HD 60364-7-710 standards.
7. The device lifetime is 10 years.
8. The instrument and the accessories shall be considered as "contaminated medical waste" when having expired its service life.
9. All the functions of the device should be checked before use.

1.1. Intended Use

Biçakcılar Luxline LED Surgical Light is used for assisting medical personnels at operating/ recovery rooms for achieving surgical operations, recoveries and diagnosis procedures by illuminating environment and patient body when the electricity goes off.

2. Description and Technical Specifications

2.1. General Specifications

Spring arms : 60° right, 60° left; 910 mm length,
50° downward, 45° upward,

1180 mm action distance

Hanging bridge : 360°

Action diameters : 910 mm

Lamp heads : 1 w/single LED head,
270° rotatable

Lamp head diameters : 710 mm.

Lighting Unit : LED

Life of lighting Unit : >50.000 hours

Maximum lighting : 160.000 lux

Adjustment Lighting Power : Gradually

Color Temperature : 4000 K° / 4500 K° / 5000 K°

Color Rendering Index (CRI) : 95

Light field diameter : 180 - 250 mm.

Adjustment of light field diameter : Gradually

Depth of Illumination (L1+L2) (20% Ec) : 1200 mm.

Mains power/ generator input : 230 V, 50/60 Hz

Automatic switching to battery circuit from mains circuit

Battery inlet voltage : DC12V - 12Ah

Battery recharge time : Between 15/ 18 hours

Battery Type : Dry type, no need for maintenance

Illumination time : 3 hours

Applied standards :

IEC60601-1, IEC60601-1-2, IEC60601-2-41,

Class I, suitable for continuous operation, transportable device.

The performance values are the values of the product powered by 230VAC 50Hz mains. If the product is powered with different values than those of the mains, deviations in performance values can be seen.


Figure 1

3. Packaging, Storage and Transportation

Surgical lights are shipped with a special Styrofoam packaging material after being placed inside cardboard box. The device consists of fragile compound, so protect the device from instantaneous shocks during storing and transporting and do not drop the packages. Packaging material is made of recyclable materials. Disposal in suitable waste containers is recommended.

Make sure that the surgical lights are well protected during transport against rain, snow, mechanical collisions and severe shakes. Do not store the surgical lights outside.

Appropriate storage conditions are as follow:

Warning!

The device shall not be transported as handled from its lower and upper pipes and the lamp head. Otherwise balance problems may be originated.

Warning!

The battery fuse should be removed for long-time storage and transport otherwise the device may be damaged. The batteries may be depleted and can not be used.

4. Assembly And Installation


4.1. Installation of Body

Move the light body to the intended area of use. Fix the body by pressing the brake pedal of the rulettes w/brake.

4.2. Assembly of Upper and Lower Pipes


Firstly, lower pipe is passed through lower pipe housing. Then, lower pipe is located on chassis by orienting pin holes. After that, tighten M10 screw with spring washer to fix lower pipe to chassis. Finally, tighten 2 M8x8 screws on lower pipe by orienting their screw housings.

For upper pipe installation; firstly, upper pipe is inserted into lower pipe. Then, tighten 3 M8x8 by orienting their screw housings to fix upper pipe to lower pipe.

4.3. Assembly of Spring Arm

Connect the spring arm to the upper pipe as seen on the figure. During this operation, pay attention that the washer on the upper pipe completely sets.

4.4. Adjusting the Spring arm

4.4.1 Adjusting the Balance

For the Lighthouse to be carried in a balanced manner by the spring arm, adjustment of balance on the spring arm must be made. This adjustment is made as follows:


Adjustment should be made with a 5 mm Allen wrench through the adjustment hole below the part where the spring arm is connected to the rotary arm.

1. If the Lighthouse is automatically lifted;

Adjust the Allen wrench by turning it in the negative (-) direction.

2. If the lighthouse automatically falls down;

Adjust the Allen wrench by turning it in the plus (+) direction.


4.4.2. Adjusting the Height

In places where the ceiling height is not sufficient, spring arm may be adjusted for the good height to avoid hitting the spring arm to the ceiling.

Height stop located on the edge of the spring arm towards the lighthouse must be adjusted using 5 mm Allen wrench.

Rotate the Allen wrench to (-) to lower the stop.


Rotate the Allen wrench to (+) to lift the stop.

4.5. Assembly of the Light Head

Apply a thin layer of silicon grease on the shaft on the tip of crescent pipe before you install the lamp head.

Remove the cam piece before you mount the lamp head on the spring arm.


Remove the bolt (3) on the retaining ring. After you remove the bolt, lift the (1) retaining ring by opening its cams (2) with a screw driver and attach the lamp head. Place the cam into its house, remove the retaining ring downward and tighten the bolt (3).

4.6. Removing and Affixing the Sterilisable Handle

If you want to remove the handle, press the locking pin located on the upper edge of the handle and pull the handle down while holding the pin down.

To replace the handgrip, attach the handgrip on the aluminum part and slide it into place from bottom to top. Make sure the locking pin snaps into place. Turn the handle slightly to make sure that the handle is fully seated.

Do not touch the resulting aluminum part after removing the handle. Do not move the head without the handle.

Sterilisable Handle

Attention!

Use the handle provided by the Bicakcilar,

Warning!

Maximum sterilization temperature of the handle is 134°C and maximum sterilization duration is 7 minutes.

Warning!

Sterilisable Handle can only be used by a sterilized personnel.

Warning!

Sterilisable Handle must be sterilized every time before use.

Warning!

Check the handle for cracks and whether the lock works properly after every sterilization and before every use and make sure that you lock the handle to the lamp ("Click") during surgery to avoid the risk of falling on the floor.

5. Use

5.1. Preparation for Use

Before using the light, ensure that the installation and positioning processes indicated in Section 4 have been completed. Control the fuses of the device (See Section 6.5). Insert the plug of the device into an earthed socket.

Warning!

The device should be recharged for 24 hours before use.

Warning!


Before using the operating lamp, make sure that the lamp is in good working order. Do not operate the lamp if there is a broken part on the lamp or if the lamp glass is broken or cracked.

Warning!

Do not direct the lamp from the hand grip or head. Otherwise, the lamp may tip over.

5.2. Description of the Control Panel

- Input:** Shows the mains power.
- Charge:** Shows the charge energy.
- Discharge:** Shows that the device is operating via battery.
- Poor Battery:** Indicates low energy in the battery. The device must be recharged.
- Empty Battery:** Indicates that battery energy is totally depleted. The device gives continuous alarm.

Warning!

The lamp does not work when the battery is totally depleted. Be sure to RECHARGE THE DEVICE.

- Excessive Heat:** Shows the heat inside the device.
- Battery Charge level:** Shows how much the battery is charged. The battery charge level LEDs show that the battery is recharged by turning on and off level by level. When operating via battery circuit (during discharge), the charge level LEDs show the charge level of the battery.
- Excessive Heat Alarm Reset:** The device gives intermittent alarm when the excessive heat alarm reset LED flashes. To turn the excessive heat alarm reset, alarm mute button must be pressed.

5.3. Adjusting the direction of light

Holding handles are used to position the light field and adjust the direction of light on the surface to be illuminated. The light can focus after the lamp is turned on.

5.4. Control keys


Membrane keypads on both sides of the head can be opened and closed and functions can be controlled.

Power status indicators

- LED, RED, flashing: Input voltage is not correct.
- LED, GREEN, flashing: Lamp usage time exceeds 20,000 hours
- LED, GREEN, steady light: lamp working.

Endo light status indicators

- LED, GREEN, steady: Endo light on.
- LED, off: Endo light off.

5.5. Opening / Closing the Lamp


To turn on the light, press the "ON / OFF" button. If the "ON / OFF" key is pressed again while the head is open, the light turns off. When the lamp is working, the GREEN LED status light on the membrane keypad illuminates.

5.6. Adjustment

The brightness of the light can be dimmed by pressing the DIM - button and increased by pressing the DIM + button. The brightness of the light can be adjusted from 100% to 20% in 10% steps. The blue LED indicating the light intensity level is indicated by the indicator.

The light intensity adjustment function can also be controlled by the non-contact sensor, which is integrated in the bottom surface of the lamp head. If you move your hand closer to this sensor, or if an object of similar size approaches this sensor, the sensor will adjust the programmed light intensity levels. The sensor can be triggered at a distance of 10 cm or less. For the change of light intensity, the sensor's detection time can be set to a period more than 500 ms.

The sequence of changing the light intensity of the sensor is shown in the following graphic.

Follow the diagram clockwise.

Start - Light intensity area 80% - Sensor - Light intensity area 60% - Sensor - Light intensity area 40% - Sensor - Light intensity area 20% - Sensor - Light intensity area 100% - Sensor

After the cycle is completed, the light intensity area reaches 80% again and the cycle continues.

5.7. Large / Small Lighting Focus

It is possible to switch between large and small focal diameters thanks to the light focus button.

5.8. Endo Light

The Endo key turns on the green light for operation with Endoscopy. This light is a light used by the user for tracking purposes. When the endo light is active, the green LED on the button lights up.

When the endo light is on, the small or large light field adjustment keys do not function. The light intensity setting and sensor are deactivated when the endo light is on.

To turn off the Endo light, the Endo key can be pressed or the lamp head can be turned off completely.

5.9. Activating / deactivating the sensor

The sensor is activated in normaly. The sensor can be deactivated using a key combination.

To disable the sensor:

- 1. Press and hold "DIM -"
- 2. Press the "ON / OFF" button
- 3. Release both keys

The lamp head will open and the sensor will become inactive. This will be maintained.

To reactivate the sensor, first turn off the lamp head. Then, follow steps 1 - 3 above. The sensor will be activated again.

5.10. Adjusting the color temperature

Press the CCT key to adjust the color temperature.

3 different color temperatures are defined (4000K / 4500K / 5000K). Each time you press the CCT key, the color temperature changes between these 3 values.

Since the hood is opened, the current color temperature value is 4500K.

The order continues as 4500 K - 5000 K - 4000 K - 4500 K - ...

6. Maintenance / Cleaning And Disinfecting

Surgical lamp is a product developed and produced with LED technology. The LEDs used are maintenance-free light sources. No wear and tear for the electronics. Maintenance is therefore limited to simple cleaning of the lamp head.

The lamp head is designed to be maintenance-free. During the life of the product, it is recommended to lubricate all rotational connection points with an acid-free lubricant every 2 years to ensure proper system operation.

Element	Period
Checking the functions and light intensity	Functions and light intensity must be checked before every use. A full functional and performance control must be carried out in every 2 years by an authorized personnel of Bıçakcılar.
Checking cable/wire	Functions and light intensity must be checked before every use. A full functional and performance control must be carried out in every 6 years by an authorized personnel of Bıçakcılar.
Full system controls	Functions and light intensity must be checked before every use. A full functional and performance control must be carried out in every 10 years by an authorized personnel of Bıçakcılar.

Attention!

Surgical lights must be maintained and/or serviced only by authorized personnel of Bıçakcılar.

Attention!

Make sure to only use the spare parts supplied by Bıçakcılar.

Warning!

Use the handle to adjust the position of the light.

Attention!

Avoid looking directly at the light. Looking directly at the light can cause excessive dazzle and eye damage..


Attention!

Dismantling some of the compounds may have negative effects on operation and safety (e.g. during maintenance of electric power source, rotating arm mechanism and spring arms).

6.1. Cleaning/Disinfecting

The surgical lamp should be cleaned regularly according to the rules of the organizer using the product. All disinfection measures are carried out by wiping the product by applying the disinfectant from a dosing unit on a damp cloth. The lamp housing must be cleaned/disinfected with a soft cloth. The surfaces to be cleaned should be wiped with a damp cloth. Wiped surfaces are air-dried. The contaminated surface should be wiped off quickly with a damp disposable cloth moistened with disinfectant. In principle, the effectiveness of the cleaning procedure used (eg according to the products and procedures included in the list of disinfectants and disinfection procedures that the Robert Koch Institute (RKI) has tested and approved by the German Society for Hygiene and Microbiology (DGHM) or Association Française de Normalisation (AFNOR), must be validated.

Light-transmitting surface of the lamp

Warning!

Do not use abrasive chemical agents to clean the operating lamp.

Do not apply the cleaner/disinfectant directly on the lamp head, be sure to apply with a soft cloth.

Attention!

Tested disinfectant: TEGO-CID
(Microbicidal amphotensides and Cation-active biocide aqueous solution)
Within 100gr: %10 didecyldimethylammonium chloride
%10 amino acetic acid

6.2. Cleaning the Transmissive Surface

To clean the light-transmissive surface of the lamp, disinfectant must be wiped with a soft microfiber cloth.

Attention!

Do not apply disinfectant to the light transmissive surface of the lamp.

Attention!

Contact with hard objects or corrosive chemicals on the lamp glass will damage the lamp glass.

Attention!

The lamp glass must be wiped with a soft microfiber cloth. Otherwise, the lamp glass will be damaged.

6.3. Handle Cleaning and Sterilization

Basic knowledge

It is absolutely necessary to effectively clean the handle. It is your responsibility to sterilize the product. It is the user's responsibility to sterilize the handle. The disinfectant suitable for disinfection and the sterilization method specified in the instruction manual must be applied within this scope. In addition, the hygiene rules of the hospital / clinic should be followed.

Cleaning / Disinfection

Sterilization

The hand grip can be sterilized by autoclave.

Recommended parameters;

For steam sterilization; 134 ° C, 2.3 bar, 18 minutes

(or for steam sterilization; 121 ° C, 1.3 bar, 25 minutes - 30 minutes)

The wear of the sterile hand grip is normal. Therefore, it must be checked regularly for material defects and wear.

Attention!


Do not use damaged hand grips the user is responsible for the hygiene and disinfection required.

Attention!

Hand grips are not sterile upon delivery. They should be sterilized with autoclave before the first use.

7. Labeling Information

8. Disposal

The surgical lamp does not contain any harmful material. At the end of the product life cycle, surgical lamp parts may be separated in accordance with the legal regulations of the country where the product is used. Please make sure the materials are separated correctly.

- Circuit boards can be recycled properly.
- Other elements must be separated and disposed of as recyclable waste

9. Technical Specifications


Intensity of illumination Ec @ 100cm	160,000 Lux + 0 / - 10%
Irradiance E (typ.)	527 ± 20 W/m²
Quotient E/Ec (typ.)	3.58 ± 0.15 mW/m²/lux
Light field diameter - large light fieldd10	380mm
Light field diameter - small light field d10	280mm
Colour temperature	4000K, 4500K, 5000K
Colour rendering index Ra	95
Illumination depth (L1 + L2)(@20% Ec)	1200mm
Remaining illumination no tube / with 1 mask	48%
Remaining illumination no tube / with 2 masks	50%
Remaining illumination with tube	100%
Remaining illumination with tube / with 1 mask	48%
Remaining illumination with tube / with 2 masks	50%
Power supply	24VDC ± -10%
Power consumption (For each title)	Max 70W
Operating life	50,000 h

10. Trouble-shooting & Service

Error symptom	Possible causes	Resolution
Light does not light up		Analyse error symptom in accordance with this table
The light flashes	Drop in voltage is too great through the supply line in the support system.	Consult your fitter or technician
	The voltage is set too low at the power supply.	Consult your fitter or technician
Status LED flashes red	Lamba voltajı izin verilen sınırların dışında	Consult your fitter or technician
Status LED flashes green	The light has reached 20,000h of operation.	Consult your fitter or technician
Sensor does not work	Sensor is deactivated	Activate sensor
	Protective surface is dirty	Clean protective surface
The light brightness changes constantly	Protective surface is dirty	Clean protective surface
	Sensor is defective	Deactivate sensor

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