



Features

MCB350 Double Crank Adjustment Manual Bed

Beds with hand crank adjustments of back rest and knee rest, ratchet mechanism adjustment of foot rest, height adjustment by hydraulic system and come with full PP(ABS) side railsoft-drop mechanism. Also PP (ABS) Head and Foot board provided with bed for hygean purpose. It is designed to respond to short and long term basic care needs and to be used for patients who needs caregiver supervision.

Suggested area of use : Home Care, Hospitals, Nursing Home

Supports Semi-fowler and Fowler positions. Back rest adjusted up to 0°-70° and knee rest 0°-35°, height adjustment 45-65 cm. For long time use, fowler/sitting position prevents patients from getting bedsores. With fowler positions backrest and kneerest to be adjusted simultaneously for comfortable positioning. Also with fowler position, patients will have better blood circulation which will provide faster recovering. Patients can experince an ideal sitting position with adjustment of backrest to 0-70°, therefore activities like eating, using bedpan, watching tv, reading book or using laptop to be done comfortably.

Side Rails: Fold-away PP side rails developed to provide comfortable getting in and exiting of the bed and safe recovery as they get the full length side covered. Its stable constuction and ergonomic holding points help patients to exit beds easily and safely. It contains 2 parts of side rail on each side, can be controlled seperately. It's developed for those having problem while sitting and getting up as well as providing full safety against risk of falling out . It moves up and down smoothly thanks to its gas spring mechanism and "soft drop" feature. While It is up ward position, make sure it's clicked and locked with its automatic lock mechanism. It includes angle indicator on side rails, you can adjust the backrest to desired position considering situation of the patient on bed. The PP side rails are washable, autoclavable, resistance to cleaning chemicals, light-weight, durable and minimized risk of injuries thanks to smooth finishes. Overall, the bed has appearance of professional and top-level use.

Head and Foot Panels: Made of polypropylene (PP) medical grade plastic for hygen and cleaning purpose. Washable, autoclavable and resistance to cleaning chemicals. They are designed as detachable for emergency intervenes as well as making up the bed. Height of the panels is 33,5 cm from mattress platform.The boards easily fit on sockets placed at four corner of the bed.

Mattress Platform: 4 sections, 3 adjustable platform consist adjusments of backrest 70°, kneerest 35° by hand crank and footrest 45° by ratchet mechanism. The platform has 85cm width and 190cm length for comfortable recovery.

Painting and Coating : Construction and frame surfaces are coated with epoxy powder paint at high degree ovens. Durability, resistance to scratching and smooth surface at high quality is under guarantee. Therefore, It keeps the look like the first day. Eco-friendly smooth surfaces help you take care of your beds at ease.

Construction (Body) : Made of steel construction. With gas metal arch welding technique, It provides strong and durable body.

Additional Option: Legs and boards of our beds are demountable, ther

efore It ensures maximum use of space during storage and transportation. There is space under beds needed for placing lifts and overbed tables.

Adjustment Control : Adjustments made via hand crank control. Beds equiped with mechanism that helps caregivers adjusting the bed smoothly and in silence. Cranks gears have self-lock system. Hand cranks can be foldable and hidden while not in use. That prevents injuries from accidental and unwanted hits.

Mobilisation: If you need to move a patient with the bed, then It is a must to equip beds with castors. Especially for hospitals the mobility of beds is the basic need, yet for home care use we provide options for the beds to be without castors. Also more cost-effective option, beds without castors can be a good choice.