





Innovation Deluxe Dialysis



The Innovation Deluxe Dialysis couch has been designed for use within specialist renal dialysis environments. The couch features auto-electric CPR and fully electric profiling with a 2-way tilt facility to ensure a smooth transition when positioning patients for dialysis procedures and treatments.

Models

MG3690

MG3490

Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large central locking wheel design with steering facility Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large independently braked wheel design

State
Clear
Head
Pape
Multi
Adjus
Chois option


Model MG3490 with large independently braked wheel design

Deluxe Features

- Electric height range from 47cm to 97cm
- Electric trendelenburg tilt to -/+15°
- Electric profiling backrest angle from 0° to +80°
- Electric footrest angle from +10° to -35°
- Safe working load and lifting capacity of 250Kg (550ll
- Auto-electric CPR function via han
- Auto-intelligent operation
- Vertical lifting design to minimise space and provide maximum stability throughout all profiling positions
- Smooth moulded base cover minimises dirt traps an facilitates easy cleaning
- Anti-bacterial powder coated frame and base cover
- Large Tente[®] Ø125mm castors
- SEERS ComfortPlus[™] Upholstery features anti-micro vinyl incorporating pressure relieving design for addi patient comfort with a foam depth of 60mm
- Standard width upholstery of 70cm
- Clearance beneath frame for mobile hoist
- Head support cushion fitted on Deluxe models
- Paper roll holder fitted on Deluxe models
- Multi-adjustable arm rests for patient comfort
- Adjustable upholstered foot board for patient comfor
- Choice of upholstery colours, for full range of colour options see page 51



Fully profiled for patient comfort



Height range 47cm to 97cm

	Optional Accessories 6054 Telescopic IV pole 6057 Battery back-up facility 6241 Patient lap tray
os)	See page 48-49 for additional information on our range of optional accessories.
bial tional	
t	8 8 8 8 8
١	Low patient transfer height



2-way trendelenburg tilt facility to +/-15°