
Desktop Pulsating Vacuum Steam Sterilizer

TM-20/24DV



Use

The table type sterilizers are safely auto-controlled, which are designed for the dental clinics, hospitals, school lab, scientific research institutions and laboratory. And mostly suitable for sterilizing of surgical, dental and eye instruments, glassware medicine, culture medium and biological dressing, food, etc.

Characteristics

- Standard Class B with three times vacuum and drying, the remaining temperature of sterilized instruments less than 0.2%.
- LCD display working status, touch keys.
- The vacuum data can reach -0.8bar, it fits various instrument, including package, unpackaged, solid, hollow, multi apertures and inset pipes. It makes sure the instrument can be sterilized completely.
- Fully computer controlled LCD display bright show. The interface is modularized panel for easy operation.
- It is installed with BOWIE&DICK which measures the penetration of water steam.
- It has vacuum test which can make sure the ability of vacuum testing.
- Fast and independent steam generator.

- Built-in mini printer which can re
- Equipped with Safety valve, Safety lock door system, Pressure lock system, Pressure or temperature overload protection, Alarm system.
- It inactivates the HIV, HBV, BSE, and bacillus efficiently.
- This product is designed for open water tank, water purification, sewage tank separation, convenient and intuitive.
- The cycle cannot be started if the door is not sealed. Once the cycle is started and running, the sealed door is locked and cannot be opened to ensure operational safety.

Technical Data

MODEL	TM-20DV	TM-24DV
Technical Data		
Sterilizing chamber volume	18L φ250×360mm	23L φ250×470mm
Design pressure	0.3MPa	
Design temperature	150°C	
Max.working pressure	0.25MPa	
Max.working temperature	138°C	
Timer range	0~99min	
Heat average	≤±1°C	
Power	AC220V.50Hz/3KW	AC220V.50Hz/3KW
Dimension (mm)	480×680×455	480×680×455
Exterior package dimension (mm)	580×790×580	580×790×580
G.W/N.W	56/47Kg	57/53Kg

Details

JIBIMED



Turnkey Supplier

Siège: 27, Old Gloucester Street, LONDON, WC1N 3AX, (UK) - Tél. +44 20 3398 1339
Bureaux: 15 rue de la Garenne 27950 Saint-Marcel (France) Tél. +33 1 87 66 35 60

<https://ibsgroupe.com> – Email: dircom@ibsgroupe.com